


# The GLEN REPORT

Newsletter

Volume 31.1 | May 2019

## Introducing the New GlenResearch.com

Since 1986, Glen Research has been providing the highest quality reagents for DNA and RNA oligonucleotide synthesis, backed by rigorous quality control and analysis. The redesigned Glen Research website is part of our commitment to better serve our customers by equipping them with the necessary tools for the entire oligonucleotide synthesis workflow.



United States | Sign In | Register | Quick Order | 0 Items (\$0.00)

Enter Document Name, Product Name, or Keywords

**dA-5'-CE Phosphoramidite (10-0001)**  
3'-Dimethoxytrityl-N-benzoyl-2'-deoxyAdenosine,5'-[(2-cyanoethyl)-(N,N-diisopropyl)]-phosphoramidite

**1. Size: Select a Pack Size**

0.25 g (\$75.00)	0.5 g (\$150.00)	1.0 g (\$300.00)
------------------	------------------	------------------

**2. Format: Select a Format (See Compatibility Chart):**

			
Expedite (E) 1oz Boston Round	Mermaid (M) 30mL Boston Round	ABI 15ml Septum	ABI 6.5ml Septum

**3. Catalog Number: Select Pack Size and Format to Reveal**

**4. Quantity: Select a Quantity**

- 1 +

[Add to Cart](#)

[How do I Request a Quote?](#)

**Products Specifications**  
Formula:  $C_{24}H_{32}N_8O_7P$   
M.W.: 857.95  
F.W.: 313.21  
CAS Number: 140712-82-9

### New Features

Continual innovation for over 30 years has generated a large catalog of monomers, supports, and synthesis reagents. Robust features have been implemented to help you find the most relevant products for your needs:

- Search for products by function, application, or synthesizer
- Filter by various options such as synthesis step, reagent type, protection type, and more

- Select the right format for your experiment using our compatibility guide
- View reagent details, technical documents, dilution & coupling data, FAQs, and related items from product pages
- Submit your customer service, technical support, or new product requests online

*Continued on Page 2*

Introducing the New GlenResearch.com 1

5-Bromo- and 5-Iodo-Pyrimidine Nucleosides; Efficient Oligo to Protein Photo-Cross-linkers (Part 1) 3

New Product - dmf-dG-5'-CPG 11

New Product - Methacrylate C6 Phosphoramidite 13

New Product - 2'-F-5-Me-U-ANA 15

**Glen Research**  
22825 Davis Drive  
Sterling, VA 20164  
Phone: 703-437-6191  
support@glenresearch.com

[glenresearch.com](http://glenresearch.com)



## Introducing the New GlenResearch.com (cont.)

Continued from Front Page

### Online Ordering

Ecommerce functionality allows for quick and easy ordering. Create an online account for additional benefits:

- Customer-specific pricing
- Fast check out
- View order history
- Quickly purchase previously ordered products
- Easily manage shipping and billing addresses

The screenshot shows the 'My Dashboard' page. On the left is a navigation menu with 'Account Dashboard' selected, and other options: 'My Orders', 'Address Book', 'Account Information', and 'Newsletter Subscriptions'. The main content area is divided into three sections: 'Account Information' with 'Contact Information' (Glen Research, Edit, Change Password) and 'Newsletters' (You subscribe to 'General Subscriptions', Edit); 'Address Book' with 'Manage Addresses' and two address entries: 'Default Billing Address' and 'Default Shipping Address', both with 'Edit Address' buttons.

### Resources & Support

We are eager to share information in ways we hope will benefit you and the industry as a whole. The Glen Research website is your gateway to our team's knowledge and expertise:

- Glen Reports
- Safety Data Sheets
- Application Guides
- Literature Highlights

We want to thank you for your continued support. Explore the site—and tell us what you think. Your feedback is very important to us as we strive to develop a website that is valuable to you.

The screenshot shows the Glen Research website's 'Guides and Procedures' page. The top navigation bar includes 'United States', 'Sign In', 'Register', 'Quick Order', and '0 Items (\$0.00)'. A search bar is present with the placeholder 'Enter Document Name, Product Name, or Keywords'. Below the navigation is a menu with 'Catalog', 'Products', 'Applications', 'Instruments', 'Resources', 'Glen Report 30.2', and 'Contact Us'. The main content area shows 'Home > Resources > Guides and Procedures' and a list of guides with PDF icons: 'TBDMS-Protected RNA Guide', 'TOM-Protected RNA Guide', 'Guide to PolyPak Purification', 'Guide to GlenPak Purification', and 'Guide to Glen GelPak Purification'.

# 5-Bromo- and 5-Iodo-Pyrimidine Nucleosides; Efficient Oligo to Protein Photo-Cross-linkers (Part 1)




**Author: Dick Keys, Ph.D.**

*This Glen Report is Part One of 2; Part 2 will be released in a subsequent report.*

## Introduction

In a recent Glen Report, we described the ultrafast photoreactive nucleosides <sup>CNV</sup>K and <sup>CNV</sup>D for DNA or RNA interstrand cross-linking. Alternatively, the photo-cross-linking reactions between modified nucleosides in DNA and RNA oligonucleotides with electron-rich amino-acid side-chains of proteins have proven broadly useful in elucidating the interaction of nucleic acids with specific binding proteins, and has seen an expanded interest in studying photo-cross-linking of nucleic acids to ribonucleoproteins. In an earlier Glen Research literature review of oligonucleotide crosslinking, the halogenated nucleosides received brief mention compared with other approaches; yet 5-halogenated pyrimidine nucleosides are especially useful in site-specific cross-linking of oligonucleotides to specific protein binding sites.

In this report we will describe improved methods for the photo-cross-linking of single- and double-stranded oligonucleotides containing 5-halogenated pyrimidines with corresponding oligo binding proteins. DNA and RNA oligonucleotides containing these modified nucleosides are readily synthesized via standard machine synthesis using phosphoramidite monomers (Figure 1). Through examples, we hope to improve recognition of the special utility of 5-bromo- and 5-iodo-pyrimidine nucleosides in synthetic DNA and RNA oligonucleotides by providing background on a) the photolysis reaction mechanisms and their dependence on UV wavelength, b) the steric impact of 5-halogenated uridine containing oligos



**Figure 1. 5-halogen-modified nucleoside CE phosphoramidites**

on binding affinity, and c) conditions used to obtain high efficiency photo-cross-linking to oligo binding protein. Hopefully data presented here can dispel a misplaced notion that low cross-linking yields are obtained when using these compounds, and that their highly specific photo-cross-linking reactions will lead to additional research utilization.

Long-time readers may recall that 5-halogen substituted deoxyuridine (dU) has been useful in solving oligonucleotide protein complex structures using X-Ray crystallography. In addition, 5-Bromo-dU (as triphosphate) incorporated into nascent transcripts has found application in cell cycle investigations where antibodies to 5-Br-dU are used for immunostaining of nuclear structures. Similarly, synthetic

oligonucleotides containing 5-Br-dU are also useful as in situ hybridization probes, providing a route towards highly specific immunostaining of mRNA.

The use of short wavelength UV<sub>254nm</sub> cross-linking (UVC) of native nucleic acids has utility in discovery research, but in some applications is limited by an associated propensity for non-specific covalent cross-links between double stranded nucleic acids and non-specific coupling to numerous nearby amino acid residues in nucleic acid-protein complexes. The concurrent use of sensitizers in photochemical cross-linking expanded the utility of photoreactions by providing higher cross-linking efficiency and has been shown to increase specificity by facilitating use of longer wavelength UV radiation (UVB),



## 5-Bromo- and 5-Iodo-Pyrimidine Nucleosides; Eff (cont.)

affording fewer non-specific nucleic acid interstrand cross-links and strand breaks. For example, in the case of psoralens as photosensitizers, interstrand photo-cross-linking and photo-reversal of crosslinks facilitate precise mapping of folded RNA substructures.

Early photo-cross-linking studies with Br-dU modified oligos used short wavelength UV<sub>254nm</sub> and results were confounded due to photo-activated side reactions (see below). However, the Tad Koch lab at the University of Colorado, Boulder demonstrated that lower energy UV irradiation (UVB or UVA) could be used to minimize protein and nucleic acid side reactions and promote high-yield photo-cross-linking of halogen-nucleoside containing oligo with protein when complexed to the target protein. In model studies of the photo-cross-linking reactions, irradiation wavelength dependence and its relationship to the formation of different products were studied, and two major routes to these differing products were elucidated.

While UV<sub>254nm</sub> results in homolytic C-Br bond cleavage and generation of a highly reactive uridyl radical, the formation

of highly specific cross-linking was proposed to occur using lower energy (308nm) excitation initially yielding a  $n,\pi^*$  singlet state. Through intersystem crossing, the molecule converts to a lower energy triplet state that can oxidize a nearby electron-rich peptide moiety such that the resulting uridyl radical ion will couple to an adjacent electron-rich peptide functional group, substituting it for the bromine (Figure 2. Jablonski diagram for photoactivation of the 5-Bromo-Uracil chromophore. Reprinted from "Photochemical coupling of 5-bromouracil (BU) to a peptide linkage. A model for BU-DNA protein photocrosslinking," T.M. Dietz, R.J. von Trebra, B.J. Swanson, and T.H. Koch, 1987, Journal of the American Chemical Society, 109 (6), 1793-1797. Copyright 1987 American Chemical Society.).

Much has been learned since that early work. In designing photo-cross-linking studies, one should note that the photocoupling reaction can be compromised by using excess laser power since this can cause secondary cross-link breakage as well as protein photodamage as was shown using phage R17 coat protein.<sup>8</sup> In addition, intrastrand


crosslinks can form between 5-Br-dU or 5-Br-dC and adjacent 5'- or 3'-dG and-dA residues.<sup>9</sup> Yields of intrastrand crosslinks are much higher when the neighboring nucleoside is a dG versus when it is a dA, and Br-dC forms these cross-links much more efficiently than 5-Br-dU. Recently 300 nm UV (UVB) was used to study single strand breaks and intrastrand crosslinks in double stranded oligo containing 5-Br-dC or 5-Br-dU.<sup>10</sup> These studies indicate that Br-dC is not the better choice when designing oligo to protein photo-cross-linking studies. In addition, minimization of undesirable strand cleavage and intrastrand cross-linking side reactions can be impacted by the position of the halogen containing nucleoside in the oligo sequence (e.g. adjacent to dG), and by the wavelength of UV used. Despite the potential limitations, high efficiency photo-cross-linking using these molecules continues to be a valuable technique.

Today, we increasingly understand that folded RNA structures in concert with linear sequence in RNA-protein complexes are involved in the control of gene expression. As a result, researchers continue to seek better understanding of the molecular details of these regulatory ribonucleoprotein complexes, including the specific location of the interacting oligo and protein partners in the complex.<sup>11</sup> In this report we highlight a few discovery applications that include DNA- and RNA-protein cross-linking.

### Cross-linking BrdU modified DNA to Protein

#### Nuclear factor NF BA1

Studying NF BA1, a nuclear factor that binds to the promoter region of human apoB gene, Kardassis prepared a synthetic DNA oligonucleotide modified with five 5-Br-dU residues replacing T and corresponding to the coding sequence (-80 to -63) of the apoB gene (Figure 3. apoB Promoter



**Figure 2. Jablonski diagram for photoactivation of the 5-Bromo-Uracil chromophore. Reprinted from "Photochemical coupling of 5-bromouracil (BU) to a peptide linkage. A model for BU-DNA protein photocrosslinking," T.M. Dietz, R.J. von Trebra, B.J. Swanson, and T.H. Koch, 1987, Journal of the American Chemical Society, 109 (6), 1793-1797. Copyright 1987 American Chemical Society.**






## 5-Bromo- and 5-Iodo-Pyrimidine Nucleosides; Eff (cont.)

corresponding to cleavage of labeled strands at each substituted position. Thus,  $T_{per}/B$  yielded T strand fragments at-5, -3, +3, +4, +6, +14, +20, +21, +22 and  $B_{per}/T$  yielded strand +1, +2, +8, +13, +15, +16, +18, +19, +24.

The eighteen singly 5-Br-dU modified Operator duplexes were then used to study duplex binding affinity to lac Repressor. The Tx modified mimics (top strand modified) had dissociation constants,  $K_d$ , in the range 4.1 to  $8.0 \times 10^{-10}$  M and the nine Bx mimics  $K_d$  were in the range 4.0 to  $8.3 \times 10^{-10}$  M. Essentially in all singly modified analogs affinity was little changed from the completely unmodified operator sequence ( $K_d = 6.2 \times 10^{-10}$  M). These results confirm that substituting bromine for a methyl group on the pyrimidine ring has little steric influence on binding affinity, in accord with expectation, since the van der Waals radius of bromine (1.95 Å) is quite similar to a methyl group (1.99 Å).

Impact of Transcription activator isopropyl-β-D-thiogalactoside (IPTG): Each of these 18 5-Br-dU operator mimics was photo-cross-linked in the presence and absence of IPTG; as a transcription activator it was thought that IPTG binding to the tetrameric protein may affect concurrent binding of the double stranded operator sequence and, as a result, influence the specificity of any cross-linking reactions. Indeed, it was observed that cross-linking specificity was profoundly impacted by the presence of IPTG. In its absence some amount of cross-linking was observed using all the duplexes with their varied 5-Br-dU substitution positions. However, with IPTG present, photo-cross-linking was much more specific, occurring essentially at predominantly 5 positions of substitution, revealing the presence of five specific points of close contact between the operator and lac Repressor (Figure 6).

Thus, singly 5-Br-dU containing DNA oligonucleotides confirmed that halogen substitution for a methyl group has a minimal impact on duplex affinity for the lac repressor binding site. In addition, these modified operators revealed the identity of the specific points of contact between the operator duplexes and the lac Repressor, and also showed the impact of transcription activator IPTG in focusing crosslinking to a small subset of sites within the binding pocket.




**Figure 6. Cross-linking Operator to lac Repressor Protein. Solid bars-specific cross-link formation; Shaded bars-protein to DNA cross-linking in presence of IPTG inducer.**

### Cross-linking 5-Iodouridine and 5-Bromouridine modified RNA to Protein

Photo-cross-linking of a short hairpin RNA fragment to bacteriophage R17 coat protein was studied and the efficiency and specificity benefits of using lower energy UV irradiation of the 5-I-U nucleoside were demonstrated. Bacteriophage R17 replicase gene translation is repressed when a small hairpin RNA binds to the bacteriophage R17 coat protein. Earlier workers used homologous synthetic RNA variants to establish the binding site size and critical RNA sequence parameters and found a 21-nucleotide fragment of R17 RNA binds to phage coat protein with similar affinity to the natural sequence.<sup>16</sup> In photo-cross-linking studies in Tad Koch's

lab,<sup>17</sup> Gott and coworkers prepared R17 RNA hairpins by *in vitro* transcription from synthetic DNA templates to contain 5-bromo-U substituted for various uridine residues.<sup>18</sup> From studying photo-cross-linking of the varied 5-Br-U containing analogs to coat protein, it was clear that bromouridine in the hairpin loop was implicated, and while some of the other analogs were capable of cross-linking, these cross-linkings were strongly dependent on the position of 5-Br-U in the fragment.

Willis and co-workers<sup>19</sup> synthesized by *in vitro* transcription (IVT) native and singly substituted 5-Br-U and 5-I-U RNA variants of the close-analog 19-nucleotide R17 RNA hairpin and studied their binding affinity and photo-cross-linking to purified coat protein, using a nitrocellulose filter method and alpha <sup>32</sup>P-C labeled R17 RNA (Figure 7).



**Figure 7. Structures of the bacteriophage R17 RNA Hairpins 1, 2 and 3**

Interestingly, the dissociation constants,  $K_d$ , for the iodo- and bromo-modified RNAs bound to coat protein revealed 3- and 5-fold *stronger* binding affinity to the coat protein, respectively, than the native sequence. Thus, Willis et al. found 5-Br-U and 5-I-U labeled R17 RNA hairpins were well-tolerated within the binding site of R17 coat protein; this strong binding occurs despite that the van der Waals radii of bromine (1.95 Å) and iodine (2.15 Å) are considerably larger than the hydrogen



at ring position 5 of uracil (1.2 Å).<sup>20</sup> It was noted that 5-Br-U cross-linking yields were low, and the researchers observed that photoactivation using a low-pressure mercury lamp UV<sub>254nm</sub> versus monochromatic UV<sub>308nm</sub> from a xenon chloride excimer laser resulted in different cross-linking results; using longer wavelength UV higher photo-cross-linking yields were obtained, with less protein damage and less RNA strand scission. High energy photo-cross-linking appears to result in C-Br bond homolysis forming a highly reactive uridyl and Br radical pair, while the lower energy irradiation generates an alternative excited state<sup>6</sup>.

Thus, using the 308 nm XeCl excimer laser, the 5-Br-U containing R17 RNA hairpin 1 (5-Br-U-R17 RNA1) achieved a maximal 40% coupling yield with R17 coat protein, while the 5-I-U-R17 RNA2 achieved 80% coupling yield in 5 minutes of irradiation.

The researcher's interest in further reducing photo-generated side products prompted the use of longer wavelength light to activate and cross-link the 5-I-U nucleoside modified R17 hairpin analog. The time course of photo-cross-linking

using 325 nm helium cadmium (HeCd) laser irradiation of 5-I-U-RNA 2 to the R17 coat protein is shown in Figure 8. SDS PAGE analysis revealed up to 94% coupling yields with negligible amounts of side products even upon extended irradiation.

To establish the identity of the cross-linked amino acid in the RNA-coat protein adduct, a 308 nm XeCl laser irradiated cross-linking reaction was carried out using the radiolabeled 5-Br-U RNA and coat protein. The mixture was ethanol precipitated and the redissolved pellet was trypsin digested. The mixture was purified by DEAE adsorption to remove free RNA and was subjected to a salt step gradient to remove peptides. Then RNA and tryptic adducts were eluted in 0.6 M NaCl, ethanol precipitated and purified by denaturing 20% PAGE followed by electroblotting onto a PVDF membrane. The RNA-peptide adduct was then sequenced via automated Edman degradation directly off the membrane, revealing that a single tyrosine residue, Tyr85, had been covalently cross-linked to the uridine residue.


In summary, this work showed that even the replacement of uridine by the larger 5-Br-U and 5-I-U in R17 RNA did not significantly perturb binding affinity for phage coat protein; surprisingly the already strong RNA-R17 coat protein binding was enhanced in the halogen modified R17 RNA hairpins. In addition, both the bromine and iodine singly halogen substituted RNAs were rapidly photo-cross-linked to coat protein with high efficiency using monochromatic laser irradiation at 308 nm or 325 nm, respectively; a single tyrosine residue in R17 coat protein was found to be cross-linked. Photo-cross-linking studies using longer wavelength UV irradiation resulted in reduced side-reactions and are more specific for 5-halogen-uridine containing RNA and DNA cross-links.

Low energy UV photo-cross-linking can be achieved using quite simple equipment. Dietz *et al.*<sup>21</sup> demonstrated similar specific cross-linking results to those achieved with the laser by excitation instead using the 313 nm emission from a high-pressure mercury lamp selected with a monochromator at 310 nm with a band pass of 20 nm. Gott and coworkers used a medium wavelength transilluminator for some experiments. Even simpler methods to obtain useful photo-cross-linking conditions are available, as described by Xue and Nicholson.<sup>22</sup>

## Photo-cross-linking in Aptamer Discovery and Characterization RNA Aptamer Development

Several labs have recognized the utility of photochemical cross-linking methods to better understand aptamer interaction with targets. An early application of photo-cross-linking to aptamer development by Jensen and coworkers used *in vitro* selection to direct the covalent attachment of high-affinity RNA ligands, aptamers containing 5-I-U, to human immunodeficiency virus type 1 Rev protein.<sup>23</sup>

A feature of the Jensen lab study was the application of SELEX methodology to include 5-iodouridine triphosphate (5-I-UTP) instead of UTP in the phage T7 RNA polymerase transcription reaction mixtures during evolution of the aptamers. The resulting high affinity RNA aptamers were used for efficient photo-cross-linking of the aptamer to HIV-1 Rev protein and, although this work used transcription of 5-iodouridine triphosphates (as opposed to synthetic oligos from phosphoramidite synthesis), the value of photo-cross-linking methodology in elucidating aptamer interaction with target protein was clearly shown.



**Figure 8. Photo-cross-linking of 5-I-U RNA to R17 coat protein with monochromatic emission at 325nm.**




## 5-Bromo- and 5-Iodo-Pyrimidine Nucleosides; Eff (cont.)

### DNA Aptamer development

The utility of directly using 5-Br-dU or 5-I-dU as a randomization monomer amidite in SELEX is complicated by the base pairing ambiguity of these two monomers. Still, there is another use for these nucleoside phosphoramidites later in aptamer studies by substituting for T residues in an already optimized aptamer. For example, Mallikaratchy et al. used SELEX to evolve aptamers from live cells.<sup>24</sup> Using Ramos cells, a Burkitt's lymphoma cell line, an aptamer, TD05, was discovered that recognizes with high affinity ( $K_d$  0.75 nM) a membrane bound heavy mu chain of IgM5 (IGHM), this latter a protein component of a B-cell receptor complex in these cells.<sup>25</sup> Then, using the fluorescent FITC-labeled TD05 aptamer, a cell binding assay was established by fluorescence activated cell sorting (FACS). Using this assay, specific binding of the aptamer to a yet unidentified target within cells could be demonstrated by competition with unlabeled TD05.

The researchers then prepared a series




**Figure 9. Modified aptamer with photoactive 5-I-dU linked to biotin via a disulfide bond. U, 5-I-dU; S, sulfur; PEG, polyethylene glycol linker**

of 5-I-dU modified TD05 aptamers (5-iodo-2'-deoxyUridine is called 5dUI in the paper) using machine DNA synthesis where 5-I-dU was substituted for thymine, including fully substituted (all T replaced) and various partially substituted TD05 analogs. Fully substituted 5-I-dU TD05 did not bind to cells, however less-substituted variants did bind strongly; an optimally labeled 5-I-dU TD05 aptamer analog was identified that contained four 5-I-dU residues (Figure 9). An analogous variant was then synthesized by solid phase synthesis that also incorporated a 3'-cleavable disulfide linked biotin tail to capture and identify the protein target.

To identify the aptamer target, the 5'-<sup>32</sup>P radiolabeled modified aptamer was photo-cross-linked to its target within cells using nanopulsed XeCl excimer laser (308 nm) irradiation. The resulting cells were lysed/homogenized, and solids were pelleted to yield crude membrane protein extract; this was detergent solubilized and separated from insoluble cell debris. The soluble fraction was incubated with streptavidin magnetic beads, washed, and the cross-linked protein mixture was released from beads by disulfide bond cleavage; mass spectrometry was used to identify four candidate protein targets among what remained a complex mixture also containing nuclear proteins. Further experiments demonstrated IGHM was a TD05 target. This effort led to subsequent studies in which the TD05 aptamer was prepared covalently coupled to light-activated photosensitizer molecule chlorin e6 (Figure 10) and, using the conjugate, these workers demonstrated specific killing of cancer model Ramos cells.<sup>26</sup>

In summary, SELEX optimized DNA and RNA aptamers have been modified with 5-iodouracil nucleoside residues substituted at one or more T or U

residues, respectively, resulting in 5-iodouracil-modified aptamers with retention of strong target binding affinity. The specific, efficient photo-cross-linking within a complex milieu of live cells has facilitated identification of a membrane bound receptor specific to these Burkitt's lymphoma cells and led to development of a cancer cell-specific aptamer-drug conjugate.



**Figure 10. Chlorin e6, Light Activated Photosensitizer**

### Cross-linked Adduct Structures

The photo-cross-linking regioselectivity of these methods as applied to many DNA and RNA-protein complexes has been used to reveal molecular details of complex formation. Today, an arsenal of halogen modified deoxyribonucleosides and ribonucleosides are activated by long-wavelength UV to yield photoadducts, as shown in Figure 11 below, written generally to include deoxyribo- and ribonucleotides. Glen Research carries many of the halogen containing phosphoramidites needed for chemical synthesis of photo-cross-linkable oligonucleotides. Because these nucleosides interact mainly with the electron rich amino acids phenylalanine, tyrosine, tryptophan and histidine, researchers may encounter nucleic acid-protein complexes in which a photoactivatable probe's protein binding site may not have a suitably located




nearby amino acid side chain with which to efficiently photo-cross-link. Even so, taking advantage of oligo design flexibility a great many examples have proven successful and instructive.

Although not the main focus of this report, it should be mentioned that thio-analogs of nucleosides have also been used as photo-cross-linking probes of nucleic acid interaction with proteins.<sup>27</sup> This includes the DNA analogs 4-thio-dT, 4-thio-dU, 6-thio-dG and the RNA analog 6-thio-G. Possible advantages of using the thio derivatives for photo-cross-linking is their similarity to the natural structures and the availability of additional protein interfaces for cross-linking due to the additional photoactive analogs thioG and thiodG. Typically, thio-analogs are photoactivated at 340 nm wavelength and so, similar to irradiation of 5-Br-uracil and 5-I-uracil containing nucleosides, will result in low non-specific side reactions.

## Conclusions

A long history of photo-cross-linking using 5-bromo and 5-iodopyrimidine nucleoside analogs in synthetic DNA and RNA continues aiding researchers to understand the intimate interactions of nucleic acids with partner proteins involved in cellular control and signaling pathways. Researchers have found new applications of these molecules in aptamer design and characterization, including development of drug candidates through cell-based SELEX methods in the discovery of membrane bound receptors for potential use in targeting cancers. The use of these efficient photo-cross-linkers seems well positioned to lead to other research insights among the growing network of ribonucleoproteins and RNP complexes.

The size change resulting from bromine and iodine substitution on position 5



of the uracil ring has a minimal impact on DNA oligo binding to its receptor protein binding site. Surprisingly, halogen substitution in RNA, despite the increased size of bromine and iodine over hydrogen, did not reduce binding affinity in the examples shown. Rapid and highly specific photo-cross-linking of 5-Br-dU/5-Br-U and 5-I-dU/5-I-U residues in single and double stranded DNA or RNA oligonucleotides in association with cognate binding proteins can be achieved using longer wavelength UV irradiation (laser or bandpass filtered) centered at 308 nm and 325 nm, respectively.

Glen Research's line of molecules includes the halogenated-deoxycytidine and-deoxyuridine phosphoramidites and bromodeoxyuridine CPG support. For RNA constructs, the bromo- and iodo-uridines are available as TBDMS

protected phosphoramidites and also as the bromo-uridine, 2'-O-methyl phosphoramidite. In general, mild room temperature deprotection is required. See use instructions for dissolution, coupling and deprotection conditions for each monomer. *In Part 2 of this series, we will describe elegant studies of the RNA interference pathway demonstrating how photo-cross-linking of 5-halogen-Uracil containing oligos have been used to unravel the mechanistic details of the multiple step, multiple component ribonucleoprotein RISC complex.*

## References:

1. Glen Report 2018, **30.2**, 1-4
2. Glen Research Literature review 2011: Oligonucleotide Cross-Linking
3. Abbreviations: 5-Br-dU, 5-bromo-2'-deoxyuridine; 5-I-dU, 5-iodo-2'-deoxyuridine; 5-Br-U, 5-bromouridine; 5-I-U, 5-iodouridine; 5-Br-dC, 5-bromo-2'-deoxycytidine; 5-I-dC, 5-iodocytidine; FITC, fluorescein isothiocyanate; PAGE,



## 5-Bromo- and 5-Iodo-Pyrimidine Nucleosides; Efficient Oligo to Protein Photo-Cross-linkers (Part 1) (cont.)

- polyacrylamide gel electrophoresis; IPTG, isopropyl- $\beta$ -D-thiogalactoside; UVA, 315-400nm wavelength; UVB, 280-315nm wavelength; UVC, 100-280nm wavelength.
- W. Hendrickson, C. Ogata, *Meth. Enzymol.*, 1997, **276**: 494-523; Walsh M.A.; G. Evans G., R. Sanishvili, I. Dementieva, A. Joachimaki, *Acta Cryst.*, 1999, **D55**: 1726-1732
  - G.F. Jirikowski, J.F. Ramalho-Ortigao, T. Lindl, and H. Seliger, *Histochemistry*, 1989, **91**: 51-53; G.F. Jirikowski, J.F.

- Ramalho-Ortigao, K.W. Kesse, F.E. Bloom, *Histochemistry*, 1990, 94:187-190
- M. Ascano, M. Hafner, P. Cekan, S. Gerstberger, and T. Tuschl, *Wiley Interdiscip Rev RNA*, 2012, **3**(2): 159-177
  - T.M. Dietz, R.J. von Trebra, B.J. Swanson, and T.H. Koch, *J. Am. Chem. Soc.*, 1987, **109**, 1793-1797
  - J.M. Gott, M.C. Willis, T.H. Koch, and O.C. Uhlenbeck, *Biochemistry*, 1991, **30** (25), pp 6290-6295
  - Y. Zeng and Y. Wang, *Nucleic Acids*

- Research, 2006, **Vol. 34**, No. 22 6521-6529
- M. Zdrozowicz, P. Wityk, B. Michalska and J. Rak, *Org. Biomol. Chem.*, 2016, **14**, 9312
  - V.S. Lelyveld, A. Björkbom, E.M. Ransey, P. Sliz and J.W. Szostak, *J. Am. Chem. Soc.*, 2015, **137**, 15378-15381
  - D. Kardassis; V.I Zannis; C. Cladaras *J. Biol. Chem.*, 1990, **265**: 21733-21740
  - K.L. Wick and K.S. Matthews, *J. Biol. Chem.*, 1991, **266**, 6106-6112; T.D. Allen, K.L. Wick, K.S. Matthews, *ibid.*, p. 6113.
  - E.E. Blatter, Y.W. Ebricht, R.H. Ebricht, *Nature*, 1992, **359**, **650**, GCN4
  - J.H. Miller and W.S. Reznikoff (eds), 1980, *The Operon*, 2nd Ed, Cold Spring Harbor Laboratory Press, Cold Spring Harbor, NY
  - A. Bernardi, P-F. Spahr, *Proc Natl Acad Sci U S A*, 1972, **69**(10): 3033-3037; M. Krug, P.L. De Haseth, and O.C. Uhlenbeck, *Biochemistry*, 1982, **21** (19), pp 4713-4720; J. Carey, P.T. Lowary, and O.C. Uhlenbeck, *Biochemistry*, 1983, **22** (20), 4723-4730
  - T.M. Dietz, R.J. von Trebra, B.J. Swanson, and T.H. Koch, *J. Am. Chem. Soc.*, 1987, **109**, 1793-1797
  - J.M. Gott, M.C. Willis, T.H. Koch, O.C. Uhlenbeck, *Biochemistry*, 1991, **30**(25):6290-5
  - M.C. Willis, B.J. Hicke, O.C. Uhlenbeck, T.R. Cech, T.H. Koch, 1993, *Science*, 1993, **262**, 1255-7
  - M.C. Willis, K.A. LeCuyer, K.M. Meisenheimer, O.C. Uhlenbeck and T.H. Koch, 1994, *Nucleic Acids Research*, **22**(23) 4947-4952
  - T.M. Dietz, R.J. von Trebra, B.J. Swanson and Tad H. Koch, *J. Am. Chem. Soc.*, 1987, **109**, 1793-1797
  - Y. Xue and W.L. Nicholson, *Applied and Environmental Microbiology*, 1996, **62**(7):2221-7
  - K.B. Jensen, B.L. Atkinson, M.C. Willis, T.H. Koch and L. Gold, 1995, *Proc. Natl. Acad. Sci. USA* 92, 12220-12224
  - Z. Tang, D. Shangguan, K. Wang, H. Shi, K. Sefah, P. Mallikaratchy, Y. Li, and W. Tan, *Anal. Chem.*, 2007, **79**, 4900-4907
  - P. Mallikaratchy, Z. Tang, S. Kwame, L. Meng, D. Shangguan, W. Tan, *Mol. & Cell. Proteomics*, 2007, 2230-38
  - P. Mallikaratchy, Z. Tang and W. Tan, *ChemMedChem.*, 2008, **3**(3): 425-428
  - T.T. Nikiforov and B.A. Connolly, 1992, *Nucleic Acids Res.*, **20**, 1209-1214

Item	Catalog No.	Pack	Price (\$)
5-Br-dC-CE Phosphoramidite	10-1080-90	100 $\mu$ mole	60.00
	10-1080-02	0.25 g	160.00
5-I-dC-CE Phosphoramidite	10-1081-90	100 $\mu$ mole	135.00
	10-1081-02	0.25 g	355.00
5-Br-dU-CE Phosphoramidite	10-1090-90	100 $\mu$ mole	60.00
	10-1090-02	0.25 g	160.00
5-I-dU-CE Phosphoramidite	10-1091-90	100 $\mu$ mole	60.00
	10-1091-02	0.25 g	160.00
5-Br-dU-CPG	20-2090-01	0.1 g	50.00
1.0 $\mu$ mole columns	20-2090-41	Pack of 4	200.00
0.2 $\mu$ mole columns	20-2090-42	Pack of 4	120.00
Br-U-CE Phosphoramidite	10-3090-95	50 $\mu$ mole	98.00
	10-3090-90	100 $\mu$ mole	195.00
	10-3090-02	0.25 g	475.00
5-I-U-CE Phosphoramidite	10-3091-95	50 $\mu$ mole	98.00
	10-3091-90	100 $\mu$ mole	195.00
	10-3091-02	0.25 g	475.00
2'-OMe-5-Br-U-CE Phosphoramidite	10-3190-90	100 $\mu$ mole	240.00
	10-3190-02	0.25 g	675.00
4-Thio-dT-CE Phosphoramidite	10-1034-95	50 $\mu$ mole	165.00
	10-1034-90	100 $\mu$ mole	295.00
	10-1034-02	0.25 g	675.00
4-Thio-dU-CE Phosphoramidite	10-1052-95	50 $\mu$ mole	165.00
	10-1052-90	100 $\mu$ mole	295.00
	10-1052-02	0.25 g	675.00
4-Thio-U-TOM-CE Phosphoramidite	10-3052-95	50 $\mu$ mole	212.50
	10-3052-90	100 $\mu$ mole	425.00
	10-3052-02	0.25 g	975.00
6-Thio-dG-CE Phosphoramidite	10-1072-95	50 $\mu$ mole	177.50
	10-1072-90	100 $\mu$ mole	355.00
	10-1072-02	0.25 g	975.00

## New Product - dmf-dG-5'-CPG




Oligonucleotide synthesis typically proceeds from the 3' to 5' direction, mostly because the phosphoramidites for this direction are straight forward to synthesize. These reagents allow the synthesis of oligonucleotides for a majority of the required applications. However, there are a few scenarios where 5' to 3' synthesis reagents are necessary (Figures 1 and 2).


One of these scenarios is the generation of 3' to 3' linkages, particularly as a 3'-cap for nuclease resistance. Natural oligonucleotides are highly susceptible to nucleases, and the addition of an inverted nucleotide at the 3'-end dramatically increases resistance toward 3'-exonucleases.<sup>1</sup> The most common method of introducing such a 3'-cap is to begin the synthesis with dT-5'-CPG (20-0302) and then carrying out the rest of the synthesis as one usually would in the 3' to 5' direction. The resulting oligonucleotide would have a 3' inverted dT, and it should be noted that such a 3'-cap would also prevent polymerase extension.

A second scenario is the synthesis of oligonucleotides containing modifications at the 3' terminus that would not be possible otherwise.<sup>2</sup> Examples of such modifications include certain dideoxy nucleotides and phosphoramidite reagents that are only available as terminal labels.

A third scenario is the synthesis of attached oligonucleotides that can be extended with polymerases. This type of synthesis arrangement is the key to single cell RNAseq methods such as Drop-seq<sup>3</sup> that have been developed in recent years. In Drop-seq, the synthesis begins with mono-sized microparticles containing linkers that are stable to deprotection conditions thus ensuring the oligo remains attached to the bead (Figure 3). Briefly, reverse direction DNA synthesis is performed to synthesize



**Figure 1. Reverse DNA Phosphoramidites**



**Figure 2. Reverse DNA CPG Supports**

relatively long oligonucleotides that contain several functional regions, one of which is a stretch of T's. In a microfluidics environment, the poly T region of the oligonucleotide attached to one individual bead will hybridize to/capture the poly A tails of mRNA from one individual cell. A reverse transcriptase can then append a cDNA copy of the mRNA to the bead bound oligo. Each DNA sequence contains 2 barcode regions, one to differentiate one bead from another, and the

other to differentiate one strand from another on the same bead. As a result, the DNA on these beads can be amplified by PCR for analysis by next generation sequencing to quantify individual mRNA sequences at the cellular level.

### Use of dmf-dG-5'-CPG

For many years now, Glen Research has been providing reverse DNA synthesis reagents to the research community. All these reagents have protecting



## New Product - dmf-dG-5'-CPG (cont.)




Figure 3. Drop-seq Sequence and Addition of cDNA

groups and attachment points between the 5' and 3' positions that have been reversed. Everything else is identical; the protection groups are similar and the synthesizer protocols are the same.

Early on, Glen Research offered ibu-dG-5'-CE Phosphoramidite and ibu-dG-5'-CPG; however, several years ago, we began offering the dmf-protected version of the phosphoramidite in place of the ibu-protected version while leaving the CPG offering unchanged. Many customers who are performing reverse direction synthesis may not require the dG-CPG, such as Drop-seq researchers, but for those who require the dG-CPG as well, they are forced to use more harsh conditions necessary for removal of the ibu protecting group. As listed in the Glen Research Deprotection Guide, dmf-dG is significantly faster to deprotect when using standard deprotection conditions (Table 1). Due to this mismatch in protecting groups, Glen Research has added dmf-dG-5'-CPG to its catalog.

### References:

1. J.F.R. Ortigao, et al., *Antisense Res. & Dev.*, 1992, **2**, 129-146.
2. C.D. Claeboe, R. Gao, and S.M. Hecht, *Nucl Acid Res*, 2003, **31**, 5685-5691.
3. E.Z. Macosko, et al., *Cell*, 2015, **161**,

Table 1. Deprotection of dG

dG Protection	Temperature	Time
ibu-dG	RT	36h
	55°C	16h
	65°C	8h
dmf-dG	RT	16h
	55°C	4h
	65°C	2h

Item	Catalog No.	Pack	Price (\$)
dA-5'-CE Phosphoramidite	10-0001-02	0.25 g	75.00
	10-0001-05	0.5 g	150.00
	10-0001-10	1.0 g	300.00
dC-5'-CE Phosphoramidite	10-0101-02	0.25 g	75.00
	10-0101-05	0.5 g	150.00
	10-0101-10	1.0 g	300.00
dmf-dG-5'-CE Phosphoramidite	10-9201-02	0.25 g	75.00
	10-9201-05	0.5 g	150.00
	10-9201-10	1.0 g	300.00
dT-5'-CE Phosphoramidite	10-0301-02	0.25 g	75.00
	10-0301-05	0.5 g	150.00
	10-0301-10	1.0 g	300.00
dA-5'-CPG	20-0002-01	0.1 g	50.00
	20-0002-10	1.0 g	375.00
1 μmole columns	20-0012-41	Pack of 4	100.00
0.2 μmole colums	20-0012-42	Pack of 4	75.00
10 μmole column (ABI)	20-0012-13	Pack of 1	225.00
15 μmole column (Expedite)	20-0012-14	Pack of 1	300.00
dC-5'-CPG	20-0102-01	0.1 g	50.00
	20-0102-10	1.0 g	375.00
1 μmole columns	20-0112-41	Pack of 4	100.00
0.2 μmole colums	20-0112-42	Pack of 4	75.00
10 μmole column (ABI)	20-0112-13	Pack of 1	225.00
15 μmole column (Expedite)	20-0112-14	Pack of 1	300.00
dmf-dG-5'-CPG	20-9202-01	0.1 g	50.00
	20-9202-10	1.0 g	375.00
1 μmole columns	20-9212-41	Pack of 4	100.00
0.2 μmole colums	20-9212-42	Pack of 4	75.00
10 μmole column (ABI)	20-9212-13	Pack of 1	225.00
15 μmole column (Expedite)	20-9212-14	Pack of 1	300.00
dT-5'-CPG	20-0302-01	0.1 g	50.00
	20-0302-10	1.0 g	375.00
1 μmole columns	20-0312-41	Pack of 4	100.00
0.2 μmole colums	20-0312-42	Pack of 4	75.00
10 μmole column (ABI)	20-0312-13	Pack of 1	225.00
15 μmole column (Expedite)	20-0312-14	Pack of 1	300.00

# New Product - Methacrylate C6 Phosphoramidite



Polyacrylamide gels have long since been essential tools in nucleic acid and peptide laboratories. The flexibility in preparing polyacrylamide gels allows unparalleled control of resolution through a broad range of sizes with a large loading capacity for the purification of nucleic acids. New applications for polyacrylamide gels continue to emerge.

The copolymerization of acrylamide and bis-acrylamide crosslinker to form polyacrylamide gels is a vinyl addition polymerization reaction typically initiated by ammonium persulfate and tetramethylethylenediamine (TEMED). Ammonium persulfate acts as the free radical source and TEMED as a free radical catalyst. The persulfate free radicals react with acrylamide-containing monomers, converting them into additional free radicals, as well as polymerizing with unreacted acrylamide monomers to form the gel matrix. Including a 5'-methacrylate-labeled oligonucleotide in the reaction mixture covalently incorporates the oligonucleotide into the gel matrix.

Methacrylate C6 Phosphoramidite (Figure 1) is a modifier that can be used to attach methacrylate to an oligonucleotide using conventional phosphoramidite chemistry and subsequently incorporated into a polymer (Figure 2). Methacrylate-modified oligonucleotides have been used to generate oligonucleotide-labeled hydrogels, microspheres, microarrays, and functionalized polymeric surfaces for use in purification, hybridization detection assays, affinity capture, biosensors, and sequencing. In the first hybridization proof-of-




Figure 1. Methacrylate C6 polymerization




Figure 2. Incorporation of a methacrylate-labelled oligonucleotide into a polymer matrix

concept experiments, methacrylate-labeled oligonucleotides were polymerized onto an acrylic-functionalized slide.<sup>1</sup> Subsequent hybridization with asymmetric fluorescently-labeled PCR products confirmed the utility of acrylamide-labeled oligonucleotides in detection assays. This initial research on methacrylate-labeled oligonucleotides also confirmed that methacrylate labels are stable to standard PCR conditions and once polymerized, form polymers with high thermal stability, high density, and low non-specific absorption.<sup>1</sup>

More recently, methacrylate-labeled oligonucleotides have been used in aptamer-based hydrogels, molecular imprinting, and single molecule inexpensive FISH (smiFISH) techniques.<sup>2-4</sup> Aptamer-based hydrogels use the selectivity of aptamers to bind their targets for detection using the physical properties of gels, as well as for electrochemical detection, and colorimetric detection. In one example, ochratoxin aptamers were incorporated into a hydrogel resulting in high selectivity for ochratoxin and a detection limit of 0.51 ppb in food products.<sup>5</sup>

Molecular imprinting uses biomolecules, such as antibodies, proteins, and glycoproteins, as templates to create molecularly imprinted polymers (MIPs).<sup>3</sup>



## New Product - Methacrylate C6 Phosphoramidite (cont.)

Biomolecules are mixed with the pre-polymer, polymerized in place, and subsequently removed to create a recognition cavity. Incorporating aptamers as the biomolecule in MIPs takes advantage of the binding affinity and selectivity of aptamers. As such, aptamer-based MIPs can enhance the sensitivity of MIPs down to femtomolar concentrations as shown for thrombin MIPs.<sup>3</sup>

In smFISH, a precursor to smiFISH, gene-specific primary fluorescent probes are used to localize individual mRNA within cells. By comparison, smiFISH probes are composed of primary and secondary probes. Primary probes are unlabeled and gene specific probes. Secondary probes are fluorescently labeled (FLAPs) and bind the FLAP sequence on the primary probes. When the secondary probe also contains a methacrylate

label, smiFISH probes can be used in Expansion Microscopy (ExM).<sup>4,6,7</sup> ExM is the *in situ* polymerization of labeled molecules followed by hydration of the matrix to induce expansion of the biological structure, retaining the 3D orientation of the original biological structure.<sup>6,7</sup> Combining smiFISH with ExM provides a 2-fold increase in signal-to-noise ratio and the ability to resolve overlapping transcripts.<sup>4</sup>

In each of these techniques, methacrylate-labeled oligonucleotides are incorporated into a gel matrix and the canonical properties of the oligonucleotides provide specificity and functionality within the application. Glen Research is pleased to offer Methacrylate C6 Phosphoramidite for the synthesis of methacrylate-labeled oligonucleotides.

### Use of Methacrylate C6 Modifier

Methacrylate C6 Modifier is compatible with standard coupling times and deprotection conditions. This modifier does not contain a DMT protecting group and cannot be purified with DMT-ON purification techniques. However, the purification step can typically be omitted since only the oligonucleotides containing the methacrylate modification will be polymerized into the gel matrix. Unlabeled oligos are easily washed away after the polymerization step is completed.

### References:

1. F.N. Rehman, et al., Nucl Acid Res, 1999, **27**, 649-655.
2. J. Liu, et al., Anal Bioanal Chem, 2012, **402**, 187-94.
3. Z. Zhang and J. Liu, Small, (In Press), 201805246.
4. N. Tsanov, et al., Nucleic Acids Res, 2016, **44**, e165.
5. R. Liu, et al., ACS Applied Materials & Interfaces, 2015, **7**, 6982-6990.
6. F. Chen, P.W. Tillberg, and E.S. Boyden, Science (New York, N.Y.), 2015, **347**, 543-548.
7. E.S. Boyden, 2016, <http://bit.ly/2LgJz3d>

Item	Catalog No.	Pack	Price (\$)
Methacrylate C6 Phosphoramidite	10-1891-90	100 $\mu$ mole	110.00
	10-1891-02	0.25 g	650.00

## New Product - 2'-F-5-Me-U-ANA

2'-F-arabinonucleic acid (2'-F-ANA) oligonucleotides are cousins to the widely used 2'-F-RNA, with the only difference being an inverted stereocenter at the 2'-carbon. Like 2'-F-RNA, 2'-F-ANA has high base pairing specificity and enhanced hydrolytic and nuclease stability. Glen Research began offering 2'-F-ANA versions of A, C, G and U in 2010 (Figure 1), and since then, this line of products has proven to be popular. Over the years, there have been numerous publications highlighting exciting research using these reagents.

With the goal of enhancing the pharmacological properties of G-quadruplexes, members of the Phan Lab conducted a systematic

investigation to understand the effects of several sugar modifications on G-quadruplex structure and stability.<sup>1</sup> One of these modifications was 2'-F-ANA. They synthesized a total of twenty singly 2'-F-ANA-modified variants for two types of G-quadruplexes, a (4+0) parallel and a (3+1) hybrid. Using a combination of <sup>1</sup>H-NMR, UV absorption and CD, they found that 2'-F-ANA substitutions in anti positions were well tolerated whereas substitutions in syn positions were destabilizing.

In another investigation, Holliger and coworkers investigated several different types of catalysts from non-DNA/RNA backbones such as 2'-F-ANA in the context of the origin of life.<sup>2</sup> All

their studies involved the use of *in vitro* selection (SELEX) and previously engineered polymerases.<sup>3</sup> Using a self-cleavage strategy, they were able to isolate RNA-cleaving enzymes from an initial random library of 10<sup>14</sup> different 2'-F-ANA oligonucleotides after 13-17 rounds of selection. The most active enzyme catalyst was able to cleave RNA in a site- and sequence-specific manner with an observed rate constant of 0.058 min<sup>-1</sup>. Like DNAzymes and ribozymes, this "FANAzyme" cleaves RNA resulting in the formation of a 5'-hydroxyl group and a 2',3'-cyclic phosphate. In addition to RNA cleavage, in separate experiments, they were also able to develop FANAzymes that performed the reverse reaction, RNA

## New Product - 2'-F-5-Me-U-ANA (cont.)



ligation. The RNA ligase was able to join 3'-imidazolylphosphoryl-RNA to the 5'-OH of another RNA strand, also in a sequence specific manner. Finally, they were able to use a similar selection methodology to also discover a FANAzyme that performed 2'-F-ANA ligation. The ligases exhibited rate constants of  $0.0002$  and  $0.038 \text{ min}^{-1}$ , respectively.

Recently, members of the Chaput Lab developed their own RNA-cleaving 2'-F-ANA enzymes as biologically more stable therapeutic candidates.<sup>4</sup> Also using SELEX techniques, they were able to isolate superior RNA-cleaving FANAzymes. The most active catalyst was able to cleave RNA in a sequence-specific manner with a maximum rate constant of  $0.2 \text{ min}^{-1}$  (Figure 2). Unlike earlier FANAzymes, this catalyst follows Michaelis-Menten kinetics. The researchers confirmed that the enzyme could be re-engineered to target almost any desired RNA sequence by changing the substrate binding domains of the catalyst.

### Use of 2'-F-5-Me-U-ANA

Although 2'-F-ANA are epimers of 2'-F-RNA, when 2'-F-ANA hybridizes with DNA, it adopts a more DNA-like B-type helix structure rather than an RNA-like A-type structure. Due to this, it should not come as a surprise that 2'-F-5-Me-U-ANA, which has the extra methyl group similar to thymidine, increases duplex stability relative to 2'-F-U-ANA. To give customers more control in the use of 2'-F-ANA, we are adding 2'-F-5-Me-U-ANA to our catalog (Figure 1). Like our other 2'-F-ANA phosphoramidites, 2'-F-5-Me-U-ANA can be dissolved in acetonitrile and coupled using a coupling time of six minutes using tetrazole as the activator. No changes to standard deprotection methods are required.




Figure 1. 2'-F-ANA Phosphoramidites




Figure 2. RNA-cleaving FANAzyme from the Chaput Lab



Item	Catalog No.	Pack	Price (\$)
2'-F-A-ANA-CE Phosphoramidite	10-3800-90	100 $\mu$ mole	150.00
	10-3800-02	0.25 g	375.00
2'-F-Ac-C-ANA-CE Phosphoramidite	10-3815-02	0.25 g	200.00
	10-3815-05	0.5 g	400.00
2'-F-G-ANA-CE Phosphoramidite	10-3820-90	100 $\mu$ mole	165.00
	10-3820-02	0.25 g	425.00
2'-F-U-ANA-CE Phosphoramidite	10-3830-02	0.25 g	125.00
	10-3830-05	0.5 g	250.00
2'-F-5-Me-U-ANA-CE Phosphoramidite	10-3850-02	0.25 g	250.00
	10-3850-05	0.5 g	500.00

#### References:

1. Z. Li, C.J. Lech, and A.T. Phan, *Nucleic Acids Research*, 2014, **42**, 4068-4079.
2. A.I. Taylor, et al., *Nature*, 2015, **518**, 427-30.
3. V.B. Pinheiro, et al., *Science*, 2012, **336**, 341-4.
4. Y. Wang, A.K. Ngor, A. Nikoomezar, and J.C. Chaput, *Nat Commun*, 2018, **9**, 5067.